F2 Monobloc





650/690 MM ALSO AVAILABLE IN SILVER



4-PART MATTRESS

The F2 Monobloc provides a safe, ergonomic, and reliable patient handling solution. The independent leg design supports the patient's weight while loading and unloading, reducing strain on medics.

- Easy to use controls and independent leg operation provide fast, safe loading and unloading from the ambulance
- Large 200 mm wheels easily roll over multiple terrain types
- Color-coded left-hand and righthand control levers are user-friendly and require minimal training
- All Aluminum construction is lightweight and durable



F2 Monobloc



All I

Features

- Includes telescoping IV pole, sidearms and foot support
- Independent leg operation
- Legs unfold automatically when unloading
- Easy to use color-coded controls
- Two loading heights: 560 mm and 600 mm
- Large 200 mm wheels
- Control-end swivel wheel locks
- Loading-end fixed wheels and wheel locks
- Shock frame with contour feature
- Four three-position telescoping carrying handles

- Adjustable six-position backrest
- Fold-down sidearms
- Low maintenance design
- Aluminum construction

Specifications

	MAX	MIN
LENGTH	89 in 2285 mm	76 in 1935 mm
WIDTH	22 in 565 mm	
HEIGHT	32 in 810 mm	11 in 280 mm
LOAD HEIGHT	24 in 600 mm	22 in 560 mm
WEIGHT	94 lbs 43 kg	
LOAD LIMIT	551 lbs 250 kg	

Certifications



EN 1789

Accessories



Instrument Tray
No. 0832158



•

Three Point Fastener No. 0573081



• 4-Part Mattress No. 031-4193



FERNOINTERNATIONAL.COM | +1.937.382.1457