Immobilization Solutions



Ferno's complete line of **immobilization solutions** provide **fast and safe spinal immobilization** with almost any backboard or scoop stretcher. Ferno immobilization solutions fit both adult and pediatric patients.





Fastrap Quick Restraint System quickly secures patients to backboards and other devices

Universal Head Immobilizer is lightweight, easy to apply and adaptable to any size patient



WizLoc Cervical Collar is adjustable for patients from age two through adult

GFERNO°

Immobilization Solutions

Fastrap[™] Quick Restraint System



Quickly secure patients to backboards and other devices with the Ferno Fastrap[™] Quick Restraint System. This 10-point, six-strap system is sequentially color-coded (black, red, yellow, green and gray) and includes nylon webbing with hook-and-loop closures at the ends of the straps for fast, easy application.

- · Size adjustable from pediatric to adult
- Easy storage and rapid application
- · Compatible with almost any backboard or scoop stretcher

Universal Head Immobilizer



The Universal Head Immobilizer is lightweight and easy to apply. The baseplate fits almost any backboard and the head supports can be inverted for application on a Scoop Stretcher.

- Head straps conform to the patient's forehead
- Large ear holes allow easy patient monitoring
- · Minimal interference with X-ray, MRI, or CT procedures
- PVC Plastic-coated, closed-cell foam head supports are easy to clean, will not absorb blood or bodily fluids, and prevent bacterial growth
- Meets FMVSS 302 flammability rating

WizLoc Cervical Collar



The WizLoc Cervical Collar has a universal design that fits a wide range of patient sizes from age two through adult. The adjustable, 3-position occipital support with keyhole lockpin and self-centering mandibular support provide on-patient sizing and adjustment for increased accuracy of application and greater patient comfort and safety.

- Minimal interference with X-ray or **MRI** procedures
- Can be decontaminated and is reusable up to 10 times
- Available in high visibility neon pink or OD green

Learn more at www.fernoems.com/immobilization



FERNOEMS.COM | 877.733.0911