



Model 71 Basket Stretcher



The **Model 71 Basket Stretcher** is durable and easy to handle, making it ideal for use in almost any rescue situation.

The Model 71 provides full-body protection to patients, even when they need to be transported across uneven, hazardous ground. The spacious basket stretcher offers room to accommodate blankets and life support equipment. It is available in several versions to accommodate your specific needs.

- **Ideal for most rescue situations** and available in split and non-metallic versions
- **High-density polyethylene shell** is supported by a permanently-attached, heavy-duty aluminum frame
- **Chemical, UV, rust, and corrosion resistant**, ensuring long life
- **Molded-in runners and handholds** make the stretcher easy to handle, even over rough terrain
- **Replaceable, closed-cell foam mattress** attached inside for patient comfort



5-Section Flotation Collar attaches to the Model 71 to provide flotation assistance for water rescue



Compatible with several models of **lifting bridles**

Model 71 Basket Stretchers

Specifications



	Model 71	Model 71-S	Model 71-M
LENGTH	86 in 2184 mm	86 in 2184 mm	86 in 2184 mm
WIDTH	24 in 610 mm	24 in 610 mm	24 in 610 mm
HEIGHT	8 in 203 mm	8 in 203 mm	9 in 229 mm
WEIGHT	23 lb 10 kg	24 lb 11 kg	28 lb 13 kg
LOAD LIMIT	600 lb 272 kg	600 lb 272 kg	600 lb 272 kg

Features

Model 71

- Four 2-piece, 5-foot-long quick-release patient restraints are simple to use and permit rescuers to easily change the position of straps to fit each patient
- Brass grommets are strategically located for quick use with Model 418-1 Adjustable Lifting Bridle

Model 71-S

- Separates in half for compact storage or backpacking to a rescue site
- Quick and easy to assemble and secured with pins
- Includes same features as Model 71

Model 71-M

- Non-metallic design allows for use in areas where electrical hazards could interfere with rescue efforts
- Aluminum frame is plastic-coated and fastened with plastic rivets
- Brass lifting grommets are omitted to ensure low conductivity
- Quick-fastening hook-and-loop closures replace metal buckles on patient restraints
- Features wooden runners for maneuverability, even over rough terrain

Accessories



Soft Eye Litter Bridle
No. 0313898



Model 418-1 Adjustable Lifting Bridle
No. 0312390



Model 471 Flotation Collar
No. 0 313943



Hard Eye Litter Bridle
No. 0313897



Model 513-4 IV Pole
No. 0085480

NOT PICTURED

Model 387-A Storage Case
No. 0311800 Lettered, Orange
No. 0374804 Lettered, White Cross, Orange