



Stat Trac[®]

Cot Fastening System



The **Stat Trac Cot Fastening System** protects your crews and patients in the event of an accident. The Stat Trac has been fully tested to meet ambulance safety requirements from the Society of Automotive Engineers (SAE) and is compatible with the POWERFlex[®]+ST Ambulance Cot and the 35X PROFlexx[®] Ambulance Cot.

- SAE J3027 compliant
- **Adjustable for side or center positioning** in an ambulance, making it accessible from bench seats or captain's chairs
- **Integrated Charging System (ICS[®])** eliminates the need to carry extra batteries and begins charging cot as soon as it's locked in the Stat Trac
- **Heavy-duty track system** constructed of aluminum



Straight in-track design eliminates the need to steer the cot when loading and unloading



A simple, single-hand release unlocks the cot from the Stat Trac

Stat Trac Cot Fastening System

Specifications



COMPATIBLE COTS/SYSTEMS	POWERFlexx+ST, 35XST
INTEGRATED CHARGING	Optional
NO STEERING COT	Yes
ACCOMMODATES LOAD ANGLES	Yes
LENGTH	185-8686 in 2184 mm 185-9191 in 2311 mm 185-9696 in 2438 mm 185-103103 in 2616 mm
HEIGHT	3 in 76 mm
WIDTH	8 in 203
WEIGHT	185-8655 lb 24.9 kg 185-9157 lb 25.9 kg 185-9661 lb 27.7 kg 185-1085 lb 29.5 kg

Options

- Integrated Charging System (ICS)

Certifications


The Stat Trac has been designed to meet international test standards to increase safety for you and your patient.

- SAE J3027
- KKK-A-1822
- NFPA-1917
- CAAS GVS v1.0

Features

- **Accommodates extreme loading** and unloading approach angles
- **Head-end and foot-end, dual-locking safety mechanism** ensures a secure and comfortable ride for the patient
- **Single-lever release** and a safety stop provide security while loading and unloading

Length Guide



* 25" from seat back to nearest edge of cot required for KKK compliance.

† 10" from edge of mattress to rear door required for KKK compliance.