



You carry a lot on your shoulders, and you need to be prepared with the right equipment to address any situation. FERNO **Trauma O₂^m Bags** provide multiple configurations and easy customization, helping you organize your trauma and O₂ supplies the way you want to. Every bag is built for years of dependable use under all conditions and is backed by a limited lifetime warranty.





Removable stiffeners and inserts allow custom configurations



- Concealed non-slip shoulder
 straps keep the bag on your
 shoulder
- **Reflective labels** provide nighttime visibility
- **Replaceable bottom** protects the bag and its contents from wet ground



Trauma O₂ Bags

Model 5110 Trauma/Air Management[™] Bag II



FERNO's Trauma/Air Management Bag II carries D or Jumbo D size cylinders and provides quick access to oxygen through the end pocket or from the main compartment.

WIDTH	LENGTH	HEIGHT
11 in 27.5 cm	25 in 63.5 cm	15 in 38.1 cm

- Adjustable storage in the main compartment
- Compatible with FERNO Mini-Bags
- Tuck-away backpack harness system provides hands-free carrying
- Available in blue, red, or green

Model 5111 Trauma/Air Management[™] III



FERNO's Trauma/Air Management Bag III enables excellent organization of supplies and can hold D or Jumbo D size cylinders. The bag provides quick access to oxygen through the end pocket or from the main compartment.

WIDTH	LENGTH	HEIGHT
18.5 in 47 cm	26 in 66 cm	12.5 in 31.8 cm

- Removable back pocket for WizLoc[™] or cervical collars
- Front pocket includes elastic loops and pouches to keep supplies organized and in place
- Compatible with FERNO Mini-Bags
- Available in blue, red, orange, or green

Model 5100 Airway Management/Oxygen[™] Bag



FERNO's Airway Management/Oxygen Bag is perfect for carrying supplies. The bag features two exterior utility pouches and interior space for D or Jumbo D size cylinders. Oxygen can be accessed quickly through an outside pocket.

WIDTH	LENGTH	HEIGHT
9.5 in 24.1 cm	23 in 58.4 cm	13 in 33 cm

- Includes two removable 5137 utility pouches, one 5131 pouch, and an O, cylinder bag with handle
- Tuck-away backpack harness
 system provides hands-free carrying
- **Compatible** with FERNO Mini-Bags
- Available in blue, red, or orange

Oxygen Carry[™] Bags



FERNO's Oxygen Carry Bags are designed to help you carry your oxygen cylinders onscene. The single-piece design with feet and a squared aluminum frame prevent the bag from rolling on the ground.

	WIDTH	LENGTH	HEIGHT
D	6.3 in 15.9 cm	22 in 55.9 cm	9 in 22.9 cm
DD	8.5 in 21.6 cm	22.5 in 57.2 cm	11 in 28 cm
Е	4.5 in 11.4 cm	33.5 in 85.1 cm	9 in 22.9 cm

- Padded for protection
- **Double-zip opening** for easy access to regulator
- The Model 5120 carries "D" size cylinders. Available in blue, red, orange, or green.
- The Model 5121 carries "E" size cylinders. Available in blue or red.
- The Model 5122 carries "Jumbo D" size cylinders. Available in blue.