



With a level patient surface at all heights and dual, powered actuators that maximize stability, the iNX IT is up to the difficult task of transporting your infant patients as safely as possible. It features a 700-lb unassisted lift capacity, protecting you from injury by effectively eliminating lifting.



- Maintains level patient surface at all load heights and during loading and unloading
- Battery-powered system effectively eliminates lifting during transport, protecting operators' backs



Customizable to carry most brands or models of incubator



Capable of carrying incubators and any associated equipment with 700-lb unassisted load capacity



Specifications*



LENGTH

63-81 in | 1594-2045 mm

OVERALL WIDTH 24 in | 597 mm

HEIGHT

14-49 in | 356-1245 mm

WEIGHT 208 lb | 94 kg

LOAD HEIGHT (MAX)

35 in | 889 mm

LOAD CAPACITY 700 lb | 318 kg

POWER

36 Volts DC

FASTENING SYSTEM INLINE

Features

CUSTOMIZABLE

- Customizable to carry most brands or models of incubator
- · Adapted to your specific transport needs

POWERFUL, SAFE PATIENT HANDLING

- Dual, powered actuators activate smoothly and maximize stability
- · Automatic high-speed extend & retract reduces load and unload time
- · Handles incubators and other heavy equipment with ease due to 700-lb unassisted lift capacity
- Patient surface is level at all heights to enhance procedural care
- No-lift design protects operators backs; lift, lower, load, and unload with the push of a button
- · Easy-to-use manual backup system provides complete operation at all times

TIGHT SPACE MANEUVERABILITY

- Telescoping load frame maximizes maneuverability in tight spaces
- 360-degree locking swivel wheels allow movement in any direction at any bed height

SIMPLE, SAFE, SECURE

- Simple one-touch operation helps you stay focused on patient care
- Programmable load height adjusts to match your ambulance floor height
- · Continuous and flashing lights allow you to see and be seen in dark and hazardous conditions
- LCD display with battery charge indicator provides system status and operational cues
- · Integrated Charging System (ICS®), available with iNLINE Fastening System, eliminates battery management
- In-fastener shut-off prevents accidental operation in the ambulance
- Memory-free battery and manual back-up system keep the iNX IT ready for every call

Accessories



Three-Stage IV Pole

No. LAIV3SCL

NOT PICTURED:

Loading-End Leg Storage Net

Operator-End Leg Storage Net

No. 0822454

Directional Wheel Lock

No. FWECDL

No. 0822455



Side Lift Handles

No. FWESLH



StoNet[™]

No. 0822453

 $[\]hbox{* Certain specifications may differ based on provider incubator.}\\$