

Which Stat Trac® is Right for You?

Here's a Guide to Help You Decide.

Learn more about the new standards at http://standards.sae.org/j3043_201407/



The Ferno Stat Trac Cot Fastening System has been tested to and is fully compliant with the Society of Automotive Engineers (SAE) Recommended Practice J3027 "Ambulance Litter Integrity, Retention, and Patient Restraint." Combine the Stat Trac with the POWERFlexx+ST® Powered Ambulance Cot or the 35XST PROFlexx® Manual X-Frame Cot to protect your crews and patients from injury.

Stat Trac is a straight-in-track design that eliminates cot drift and the need to steer the cot in, minimizing operator strain on angled surfaces. A single-hand release and wide-angle entrance guide provide quick and easy cot loading. Stat Trac can accept the POWERFlexx+ST and the 35XST PROFlexx cots.

*25" from seat back to nearest edge of cot required for KKK compliance.

†10" from edge of mattress to rear door required for KKK compliance.



Model 185-LT Long (103")
Typically used in Type I AD



Model 185-MT Medium (96")
Typically used in Type I or III



Model 185 Short (86")
Typically used in Mini-Mods or Type II