



## Evacuation Stair Chairs & Seats



FERNO offers a range of **Evacuation Stair Chairs and Seats** that enable fast and safe transport of a seated passenger down stairs and over flat surfaces in evacuation scenarios.

With solutions ranging from the highly-maneuverable Model 40 to the ultra-thin Model 48, FERNO has an evacuation chair that will suit your needs.

- **Quick to deploy** in evacuation scenarios
- **Ideal for use** in high-rise office buildings or industrial settings
- **Four lifting handles** with non-slip handgrips enable operators to move the passenger safely and easily



**Excellent maneuverability** for use in confined spaces



**Wide variety of options** designed for fast and safe evacuation



**Fold** for easy, low-profile storage

# Evacuation Stair Chairs & Seats

## Specifications



	Model 40-OS Stair Chair	Model 40 Stair Chair	Model 42 Stair Chair	Model 48 Stair Chair	Model 44 Rescue Seat
ITEM NUMBER	Orange: 0714841	Orange: 0714865 Burgundy: 0714864	Orange: 0714867 Burgundy: 0714868	Orange: 0714800	Orange: 0313030 Burgundy: 0313070
HEIGHT	36 in   91 cm	36 in   91 cm	37 in   95 cm	37 in   95 cm	37 in   95 cm
WIDTH	20 in   51 cm	20 in   51 cm	19 in   49 cm	16 in   40 cm	19 in   48.25 cm
DEPTH	37 in   93 cm	46 in   115 cm	34 in   86 cm	34 in   86 cm	19 in   48.25 cm
DEPTH (FOLDED)	7 in   18 cm	7 in   18 cm	6 in   16 cm	6 in   16 cm	—
WEIGHT	14 lb   6.4 kg	19 lb   8.6 kg	19 lb   9 kg	18 lb   8 kg	2 lb   1 kg
LOAD LIMIT	350 lb   159 kg	350 lb   159 kg	350 lb   159 kg	350 lb   159 kg	350 lb   159 kg

## Features

### Model 40 Stair Chair

- Easy to maneuver in tight spaces
- 350-lb weight capacity
- Ideal for use in high-rise office buildings or industrial settings
- Stationary 5-inch rear wheels enable quick transport over most surfaces
- Telescoping handles allow front-end operators to move without interfering with the passenger's feet
- Chair lock automatically engages when chair is unfolded
- Includes chest, thigh, and leg restraints

### Model 40-OS Stair Chair

- Easy to maneuver in tight spaces
- 350-lb weight capacity
- Ideal for use in high-rise buildings or industrial settings
- Stationary 5-inch rear wheels enable quick transport over most surfaces
- Similar to Model 40, but without telescoping front handles
- Chair lock automatically engages when chair is unfolded
- Folds for easy storage
- Includes chest, thigh, and leg restraints

### Model 42 Stair Chair

- Designed to be stored in small spaces
- 350-lb weight capacity
- Functions well in confined areas such as narrow hallways and multiple-landing stairways
- Seat is made of soft, easy-to-clean vinyl that resists stains
- Two 4-inch rear stationary wheel and two 2-inch front guide wheels facilitate transport on most surfaces
- Four lifting handles with non-slip handgrips help you move passengers safely and easily
- Locking safety latch maintains chair position during use and is manually released when chair is folded
- Includes restraints for chest, thigh, and lower legs

### Model 48 Stair Chair

- Identical to Model 42, but three inches narrower
- Ideal for use in airplanes, boats, restaurants, and other areas where access is severely limited

### Model 44 Rescue Seat

- Assists in moving passengers over short distances
- 350-lb weight capacity
- Four handles aid in transport
- Can remain under passenger to be used for transfers from cot to hospital bed
- Durable reinforced vinyl construction
- Includes storage pouch that can attach to a cot frame
- Easy to clean

## Accessories



### Model 559 Wall Hooks

(compatible with stair chair models 40, 40-OS, 42, and 48)  
No. 0815590



### EZ-Glide Wall Mount

(compatible with stair chair models 40-OS, 42, and 48)  
No. 0822096



### Model 59T EZ-Glide Wall Bracket

(compatible with stair chair models 40-OS, 42, and 48)  
No. 0822073



### 385-F Carrying Case

(compatible with stair chair models 40-OS and 40)  
Orange Lettered: 0311760  
Maroon Lettered: 0311762