



EZGlide® Stair Chair



The EZ Glide® Stair Chair glides patient loads up to 500 lbs down stairs without carrying or lifting. The chair and track system carry the weight, not the operator, reducing fatigue and the risk of injury.



One Inch Track-To-Ground Clearance provides smooth rolling over carpet and rough surfaces.

- **500 lb Weight Capacity** accommodates large patients.
- **Five-Position Lift Bar and Handles** reduce the strain of lifting and provide maximum control for operators of various heights.
- **Smooth 4" Front Swivel Wheels** and extra tall 6" rear wheels provide maneuverability in confined spaces such as staircase landings.

Specifications



| | | EZ Glide |
|------------------------------------|--|---|
| HEIGHT | | 37.5" (95 cm) |
| WIDTH | Overall Seat | 20.3" (52 cm) 16.5" (42 cm) |
| DEPTH | Max Open Min Open Seat Folded | 51" (130 cm) 28.5" (72 cm) 16.5" (42 cm) 10" (25 cm) |
| REAR LIFT BAR POSITIONS | | 37.5"/44"/50.5"/57"/63.5" (95/112/128/145/161 cm) |
| FRONT LIFT HANDLE POSITIONS | | 5.4"/8.8"/12.2"/15.6"/19" (14/22/31/40/48 cm) |
| WEIGHT | | 37 lb (17 kg) |
| LOAD LIMIT | | 500 lb (227 kg) |

Features

Positive Locking Hinge Mechanism

- ensures chair is locked in place whether deployed for use or folded for carrying and storage.

Large Patient Seating Surface

- with removable ABS plastic seat panels fit various size patients.

Extending Foot-Rest

- and lower track angle increase patient's sense of security.

All Metal Frame Construction

- provides durability and long life.

Upgrade Your Existing EZ Glide

- with the PowerTraxx® conversion kit.

Accessories

STORAGE



- Secure Lock +
No. 600251006



- Storage Cover
No. 0314023



- EZ Glide Wall Mount
No. 0822096


PATIENT HANDLING ACCESSORIES



- EZ Glide Headrest
No. 0822430



- EZ Glide IV Pole
No. 0821976



- Oxy-Clip O2 Bottle Holder
No. 0821977



- EZ Glide Rear Handles
No. 0821975



- EZ Glide Headrest Pouch (pouch only)
No. 0314117



- PowerTraxx Battery Charger
No. 0832199
Battery No. 0832198

Learn more at www.fernoems.com/ezglide