



Scoop EXL™ Stretcher



The **Scoop EXL Stretcher** provides superior comfort and safety to patients in need of spinal immobilization. Two hinged, interlocking pieces allow operators to bring the two halves of the Scoop together beneath the patient to gently “scoop” them up. This design eliminates the need for log-roll maneuvers, significantly decreasing movement to the patient’s cervical spine and reducing the chance of further injury.

- **Interlocking halves** allow operators to “scoop” patients up, eliminating log-roll maneuvers
- **Head section is recessed** for proper cervical alignment
- **Radiolucent**



Double-locking aluminum hinge at both ends allows either end or both ends to come apart



Telescopic tubes allow length to be adjusted to accommodate patients of various heights

Scoop EXL Stretcher

Specifications



LENGTH (MAX)	79 in 2007 mm
INTERMEDIATE 1 LENGTH	74 in 1890 mm
INTERMEDIATE 2 LENGTH	70 in 1770 mm
LENGTH (MIN)	65 in 1650 mm
FOLDED LENGTH	47 in 1193 mm
WIDTH	17 in 432 mm
DEPTH (OPEN)	3 in 76 mm
DEPTH (FOLDED)	3 in 76 mm
WEIGHT	18 lb 8 kg
LOAD LIMIT	500 lb 227 kg

Features

- **Twin Safety Lock® mechanism** provides smooth unlocking and locking
- **Ergonomic handles** support correct lifting posture
- **Thermally-treated polymer construction** does not get too hot or cold, keeping the patient more comfortable
- **Easy-to-clean surface** is impervious to fluids
- **Narrow foot-end frame** improves handling in confined areas
- **Optional speed-clip pins** allow fast attachment or rearrangement of restraints
- **Easy to store**

Accessories



TSL EXTENDER™ Kit
No. 0832274



Universal Head Immobilizer
No. 0313855



Model 430 Restraints
Contact your sales representative for product details



Wizloc™ Cervical Collar
No. 0822074 (3 pack)



365-E Head Immobilizer
No. 0314044



Model 430P Restraints
Contact your sales representative for product details



Fastrap™ Restraint System
No. 0313916



HeadHugger™ Disposable Head Immobilizer
No. 0832092 (1 pack)
No. 0832093 (10 pack)



Model 439 Restraints
No. 0313083, Black
No. 0312660, Orange



Model 436 Speed Clip Restraints
No. 0314118 (set of 4, 5' metal buckles)
No. 0314118 (set of 4, 7' metal buckles)
No. 0314118 (set of 4, 5' plastic buckles)
No. 0314118 (set of 4, 7' plastic buckles)