

FERNO ITU Incubator Transporter Unit

Our most precious cargo needs the ultimate transfer security, CUSTOM MADE* to your specifications to carry any incubator and configuration of equipment



FOR MORE DETAILED
INFORMATION SEE OVER

Specifications

Standard ITU	Length	Width	Fixed Height	Weight	SWL
	1830 mm	545 mm	615 mm	48 kg	140kg/22st/308lb

*Bespoke engineering requirements will affect specifications and test certification

FERNO – is a world leading innovator and manufacturer of patient transport solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work.

Ferno's industry-leading ITU typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999
email sales.uk@ferno.com
www ferno.com



FERNO ITU Incubator Transporter Unit

A Specialist Incubator Transfer Unit made to your equipment specifications

The Ferno Fixed-Height ITU Trolley has been specially designed to enable operators to move incubator units around and between hospitals quickly and easily.

The trolley is suitable for inter-hospital transfer by Ambulances equipped with the Ferno Lock & Load two part locking device.

So no matter what you want or how you want it, if you need to move a neonate for specialist treatment, you can rest assured you are doing the very best for their well-being and comfort during transportation.

Incubator Transfer Unit:

- **Versatility** – the trolley is designed to carry major incubator units with space saving drop-down handles
- **Smoothly manoeuvrable** – Lightweight, high quality construction enables the ITU to be operated easily
- **Static shock shielded** – anti-static wheel protection
- **Easy to clean and maintain** – manufactured from steel with powder-coat finish the paint treatments are anti-bacterial which actively inhibits bacterial growth
- **Enhanced safety** – standard ITUs are fully loaded simulation tested and certified to BS EN 1789
- **Locked down during transfer** – compatible with the Ferno Two Part Locking System

Ask about a FERNO comprehensive service contract to further protect the performance of your equipment.



Part Nos:

005103300

ITU Standard

FERNO - ITU FPC004 09.24

call +44 (0) 1274 851 999

email sales.uk@ferno.com

www ferno.com



Ferno (UK) Limited, Ferno House, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ, England.