



The New BabyPod 2 (2021 spec) is the successor to the hugely popular BabyPod II infant transportation system, but now incorporating a number of improvements and new features previoulsy introduced on BabyPod 20.

Utilising the same lid side-opening method as the original BabyPod II, but now further strengthened for improved durability, the new BabyPod 2 also includes features such as CO2 evacuation port, improved fixing of the inner foam padding, an improved finish internally and externally, and aircraft-rated 9kN trolley fixing straps for both air and road use.

BabyPod 2 shell has been designed and dimensioned to ensure it fits inside cupboards and storage areas that BabyPod II fitted into, ensuring backwards compatibility to existing vehicles.









- Side-opening lids a split-lid opens up each side of the unit
- Easy visibility and access to the baby – small opening window flaps in the upper lids allows easy access to the baby during transport whilst maintaining temperature inside the Pod
- CO2 evacuation improves baby's air supply quality

 an external inlet point on the head-end of the Pod allows air or oxygen mix to be fed into the unit, thereby diminishing any build-up of CO2

- Easy access for clinical monitoring – large access areas for tubing & cables now allow more medical equipment to be simultane-ously used
- Revised trolley strapping system – uses high-tensile webbing and metal quick-connect buckles all tested to 9kN, suitable for both air and road transport. Straps are now attached to and accessible from the outside of the Pod, making cleaning quicker and easier
- Safer snug cocoon area for precious cargo – vacuum mattress has slots in edges to allow the internal positioning straps to be fed through the mattress, thereby allowing the straps to better secure the mattress around the baby.



BABYPOD 2 (2021 spec) Infant Transport System™



Features & Benefits

- Simple to transport easily lifted in and out of vehicles, optional carry bag also available
- Unrivalled stability using vacuum mattress and patient-positioning safety straps
- Exceptionally versatile ideal for use in hospital transfers, relief operations, emergency midwifery call-outs, neonatal revival cases, in ambulance, car or air transfers
- **Convenient to use** needs no electricity, maintain temperature by use of TransWarmer heat pads

- **Proven performance** crash tested to EN 1789 and CEN compliant. Fixing restraints tested to 9kN, for use in aircraft.
- Highly compatible will fit on all models of ambulance trolley, can be used with all major ventilators enabling the addition of improved O2 environment
- **Cost-Saving** at a cost that is less than 20% of a standard transport incubator, with reduced oxygen consumption due to its compact interior and lighter due to its carbon fibre construction

Specifications

LENGTH	1000 mm
WIDTH	420 mm
HEIGHT	280 mm
WEIGHT incl. vacuum mattress	8.5 kg
CAPACITY	2 - 8 kg / 4 - 17 lbs

Made of Carbon fibre Four Layers of Protection Carbon Fibre (F1) Shell Shock-absorbing laminated foam interior Unique design Vacuum positioning mattress Vibration absorbing Gel warming mattress

Approvals and compliance: CE and FDA approved CEN compliant, dynamically crash tested to BS EN 1789 ISO 9001:2000 Quality Management ISO 13485:2016 Design + Manufacturing of Medical Devices

Accessories



Replacement Vacuum Mattress BPOD2-002

Additional Positioning Straps (single use) BP0D2-003



Replacment Shock Absorbent Foam Set BP2-005

Additional 9kN Trolley Fixing Straps BP20-006



• Replacement Lid Set BP2-007

 Protective Storage Bag BP2-013



Part Numbers

 BP2-001
 BabyPod 2 (2021 spec), includes Vacuum Mattress, Shock Absorbent Foam Set, Infant Positioning Straps (single use), Stretcher Fixing Straps, Protective Dust Cover and 1 TransWarmer

BP2-013 Protective Storage/Caccy Bag for BabyPod 2

BPOD2-020 Hand Vacuum Pump for BabyPod vacuum mattress (BabyPod II, BabyPod 2 & BabyPod 20)

- BPOD2-002 Replacement Vacuum Mattress (BabyPod II, BabyPod 2 & BabyPod 20)
- BPOD2-003 Additional single-use Infant Positioning Strap Set (BabyPod II, BabyPod 2 & BabyPod 20)
- BP2-005 Replacement Shock Absorbent Inner Foam Set for BabyPod 2 only
- BP20-006 Replacement Trolley Fixing Strap Set, metal buckles, tested to 9kN (BabyPod 2 & BabyPod 20)
- BP2-012 Replacement protective dust cover for BabyPod 2 only
- BP2-007 Replacement Lid for BabyPod 2 only
- **TRANS06**Additional TransWarmer heating pads (pack of 6)

