

FERNO PEGASUS Trolley

Ferno Pegasus is a proven versatile ambulance trolley, designed for use on ramped or tail-lift vehicles. It remains the lightest lift-assist trolley on the market today, encompassing an effortless hydraulic raise and lower system.



FOR MORE DETAILED
INFORMATION SEE OVER

FERNO – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work.

The Pegasus Trolley typifies our passion for top quality, highly functional design. . . patient transfer products that don't just do the job, they do it **BETTER**.

Specifications

| Length | Width | Height | Weight | Load Limit |
|---------------------------------------|--------|------------------|--------|-------------------------------------|
| 1950 mm | 570 mm | 880 mm (Raised) | 65 kg | 200 kg/31st/434lbs |
| 1650 mm (With folded head end) | | 460 mm (Lowered) | | 300kg/47st/658lbs (Lowest Position) |
| 1450 mm (With folded head & foot-end) | | | | |

call +44 (0) 1274 851 999

email sales@ferno.co.uk

www ferno.co.uk

 **FERNO**

FERNO PEGASUS Trolley

Quality, durability and strength... combined with unbeatable ease of use



By no means does every ambulance trolley combine optimum quality, durability and reliability and still fewer are equally suited for use on either ramped or tail-lift vehicles.

Ferno Pegasus accomplishes all of this with ease as it has to be one of the most versatile trolleys available offering exceptional simplicity of use.

- **Smoothly manoeuvrable** – extendable push/pull handles, optional folding push-poles and four pedal operated wheel/castor locks
- **Protects the patient** – Combining the use of a pressure-reducing mattress, a comprehensive harness and strap system, and moulded trolley cotsides
- **Easy to use with PTS vehicles** – standard drop-down foot end; fully compatible with Ferno lock & load system
- **Access using four-occupant size lifts** – both the head-end and foot-end fold to shorten the trolley to be able to fit into confined spaces
- **Prevent static shock** – anti-static wheels protect against electro-magnetic shock
- **Easy to clean and maintain** – removable headrest, VELCRO®-free fully welded mattress and tracker system to aid maintenance schedules
- **Versatility** – the trolley is supplied as standard with all the main bed features that are required to treat a critical patient
- **Includes easy-use two-part mechanical lock** – supplied with floor fixings
- **Enhanced safety** – tested and certified to EN 1789, CEN compliant



Pegasus with optional folding push-poles at foot-end

Part Nos:

| | |
|-----------------|--|
| PEGBIO-STD-AB | Pegasus c/w AB-PR mattress, Biosafe restraints, Std 2-Part Lock |
| PEGBIO-LHS-AB | Pegasus c/w AB-PR mattress, Biosafe restraints, LHS-operated 2-Part Lock |
| PEGBIO-NMI-AB | Pegasus c/w AB-PR mattress, Biosafe restraints, NMI 2-part Track Lock |
| PEGBIO-SNMI-AB | Pegasus c/w AB-PR mattress, Biosafe restraints, LHS-operated NMI 2-part Track Lock |
| PEGBIO-UWN-AB | Pegasus c/w AB-PR mattress, Biosafe restraints, Unwin 2-Part Track Lock |
| PEGBIO-LHUWN-AB | Pegasus c/w AB-PR mattress, Biosafe restraints, LHS-operated Unwin 2-Part Track Lock |
| PEGBIO-SL150-AB | Pegasus c/w AB-PR mattress, Biosafe restraints, SlideLock (150 mm distance) |
| PEGBIO-SL255-AB | Pegasus c/w AB-PR mattress, Biosafe restraints, SlideLock (255 mm distance) |

Also available with webbing restraints instead of Biosafe

Folding push-poles at foot end also available as an in-build option. Request for details.

SlideLock also available with bespoke distance between lock positions according to vehicle size/requirements.

Options and Accessories:

| | |
|-----------|---|
| 009001650 | Fixed, folding IV Pole |
| 009001100 | Detachable IV Pole |
| 008551000 | Oxyclip cylinder holder for CD cylinder |
| 009003410 | Base mounted dual-strap cylinder holder |
| 083195000 | PacRac Model 274-5 |
| 082-2097 | Patient shield for trolley |
| 031346990 | Trolley dust cover, grey AB material |

FERNO - PEGASUS FPC009 07.21

call +44 (0) 1274 851 999

email sales@ferno.co.uk

www ferno.co.uk



Ferno (UK) Limited, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ, England