

FERNO EZ-GLIDE PowerTraxx™

The track chair that safely ascends and descends stairs under power – at the touch of a button



FOR MORE DETAILED
INFORMATION SEE OVER

Specifications

Height	Height rear handles extended	Deployed Depth Track open, Handles extended	Seat Width	Chair Width	Depth Folded	Weight	Power Lithium Ion Battery	Max Load Capacity
953 mm	1613 mm	1295 mm	419 mm	516 mm 610 mm (with optional wide seat attachment)	254 mm	26 kg	28V 3A	227 kg/35.7st/500lbs

FERNO – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. **EZ-GLIDE** PowerTraxx™ typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



FERNO EZ-GLIDE PowerTraxx™

Safely negotiate the trickiest stairways at the touch of a button

Moving patients safely up and down stairs can be distressing to the patient and carry substantial manual handling injury risks. Ferno has the product to end these problems – the **EZ-GLIDE PowerTraxx™**.

EZ-GLIDE PowerTraxx™ is so much more flexible than any other powered tracked chair. Capable of moving down and upstairs under controlled power, you can slow it down... turn it around... bring it to a halt... all at the touch of a button.

It's not just for stairs, either – **EZ-GLIDE PowerTraxx™** effortlessly moves over all kinds of difficult terrain, including curbs. It can even be used to move heavy equipment, as well as people.

All the weight of the chair and the occupant is borne by the track and motor... not the operators – so the risk of manual handling injuries is minimised.

- **Effortless patient handling** – adjustable speed control for safe, easy use; touch-button wireless control panel for convenient operator use
- **Flexible and versatile** – climbs and descends stairs, goes around corners, tackles difficult terrain... entirely without strain
- **High load capacity** – will transport patients and equipment loads of up to 227 kgs / 35 stone
- **Excellent patient safety** – an automatic braking system combined with options such as adjustable cushioned headrest/shoulder support and extended footrest, ensures patient safety and comfort
- **Convenient to use** – detachable battery stays charged over 20 flights of steps, just 60 mins recharge time



Without extended-width seat attachment



With extended-width seat attachment



EZ-GLIDE PowerTraxx* configurations available:

- 60-0176-011 EZ-Glide PowerTraxx (Red, Rev 2 flat belt), with rear folding handles, no IV pole - Black webbing restraints
- 60-0176-011/BIO EZ-Glide PowerTraxx (Red, Rev 2 flat belt), with rear folding handles, no IV pole - Biosafe restraints
- 60-0176-011/FP/BIO EZ-Glide PowerTraxx (Red, Rev 2 flat belt), with rear folding handles, no IV pole - Biosafe & Extended footplate

Options and Accessories:

- 75-0710-001 Additional battery*
- 75-0711-002 Spare battery charger
- *All chairs come complete with 1 battery and 1 battery charger
- 082197600 IV Pole
- 082197700 CD-size cylinder holder
- 038387100 4-point restraint harness kit
- FWE2110 Cushioned headrest/shoulder support
- 082207300 Wall bracket
- 60-0251-006 SecureLock+ in-vehicle storage bracket (crash-tested)
- 009009200 Extended-width seat attachment (extends seat from 420 to 610 mm width)
- 009009200 In-vehicle storage bracket for extended-width seat attachment
- 031402300 Chair dust cover



FERNO - EZ-Glide PowerTraxx FPC006 05.19

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



Ferno (UK) Limited, Ferno House, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.