

COMPACT POWERTRAXX

Affordable Compact Powered Track Chair for everyday use



Compact PowerTraxx Chair - bringing POWER to day-to-day patient transfers

One of the most popular Ambulance Chairs used in the UK for many years is the Ferno Compact 2 Track manual track chair, widely used every day in many situations. It became the "go-to" solution to get a patient downstairs and into the back of the Ambulance. Lightweight and compact in size, it is extremely manoeuvrable around tight bends and whilst navigating narrow stairways, reliably used both in Frontline emergency service patient transfers and PTS planned transfers and patient repatriations.

But a new standard in patient transfer chair has now been established, which is based on the Compact 2 Track chair design, but powered.

Now you have the power to eliminate lifting in all situations with the Ferno Compact PowerTraxx Powered Track chair.



Compact PowerTraxx - a robust, lightweight and compact Powered Track Chair for transferring patients up and down stairways, steps, curbs and other obstacles, thereby eliminating the need to lift, also giving power-assistance to crews when moving up ramps and inclines, and contributing to an overall reduction in overall AMBULANCE vehicle and equipment weight.

Features

- Can use as a standard carry chair as well as a Powered Track chair
- Attachable Powered Track module
- Control box connects wirelessly to motor module
- Variable and fixed-speed controls
- Feature allows chair to be moved downstairs even if run out of power
- Light to illuminate dimly-lit stairways
- Adjustable-height rear handle
- Front operator grip handles
- Removable wipe-clean seat cover
- Patient headpad
- Large size wheels
- Biosafe chest & ankle patient safety straps

With many Benefits

- Not only **POWER UP and DOWN** stairways, you can now **LOAD and UNLOAD** the patient into a vehicle using the Ambulance Ramp or rear/side steps under **POWER**, and also use the **POWER** to assist on inclines, curbs, doorsteps, and other such obstacles
- The **COMPACT PowerTraxx** chair can be used on its own as a traditional carry chair, and simply attach the PowerTraxx module as and when required to use as a powered track chair, eliminating the need to lift in most situations
- Only 23.8 kg, working hard to reduce the overall equipment weight within the AMBULANCE
- Easy-use touch **CONTROL** technology, with **PROVEN LITHIUM-ION** battery power, how you want it and whenever needed.
- Operator training minimised where crews have been trained on and familiar with the Compact 2 Track chair
- Same controls and similar in operation to the Venice PowerTraxx making transition to Compact PowerTraxx quick and easy
- Proven design, exceptionally manoeuvrable in tight spaces and difficult stairways, making the Compact chair range the most used Ambulance chairs in the UK, now available in a **POWERED** version.....



Specifications

COMPACT PowerTraxx	
CHAIR HEIGHT FOLDED (without track attached)	850 mm
CHAIR DEPTH FOLDED (without track)	290 mm
CHAIR HEIGHT UNFOLDED	1102 mm
CHAIR HEIGHT WITH REAR HANDLE FULLY EXTENDED	1485 mm
SEAT HEIGHT	485 mm
CHAIR WIDTH	450 mm
TRACK ASSEMBLY LENGTH	720 mm
TRACK ASSEMBLY WIDTH	340 mm
DEPTH OF CHAIR WITH TRACK DEPLOYED	900 mm
WEIGHT OF CHAIR FRAME (NO TRACK ATTACHED)	13.4 kg
WEIGHT OF POWERED TRACK MODULE (INCL. BATTERY)	10.0 kg
TOTAL WEIGHT OF CHAIR	23.4 kg
SWL / Max load capacity (31.5 stone)	200 kg
Compact PowerTraxx in-vehicle storage bracket weight	3.9 kg

Standards:

- Compliance to EN 1865:2015
- Compliance to Medical Device Regulation (EU) MDR 17/745

Part Numbers

071481801 Compact PowerTraxx chair c/w 1 battery

Accessories



75-0710-001 Spare 5A Milwaukee M28 Battery

75-0711-002 240v Milwaukee battery charger

60-0060-001 12v Milwaukee battery charger



2-Part In-Vehicle Storage Bracket

009103900