
Product Catalogue 2012







Product Catalogue 2012

BURN DRESSINGS

Our Water-Jel Burn Dressings stop the burning process, reduce burn progression, cool the skin, protect against airborne contamination, and help relieve pain all in one simple step, with no special training and no pre-cooling required. These sterile dressings are soaked in a water-based/water-soluble, bacteriostatic and biodegradable gel and are nonadherent to burn injuries. They are composed of a medical grade, nonwoven polyester carrier. All of our dressings can be used on superficial, partial and full thickness burns. Water-Jel Burn Dressings are available in six convenient sizes and have a 5-year shelf life.



FIRE BLANKETS

Water-Jel Fire Blankets offer one of the quickest, easiest and most direct methods of administering emergency first aid burn care. When placed on a victim, our gel-soaked blankets cool the skin, protect against airborne contamination, and help relieve pain all in one simple step, with no special training required. They can also be used to extinguish flames on a victim and smother small fires.

Our Water-Jel Fire Blankets are composed of a 100% worsted-wool carrier with an intercellular weaver capable of absorbing up to 13x its own weight. The wool is soaked in a water-based, water-soluble gel that is bacteriostatic and biodegradable. The Water-Jel Fire Blankets are available in three sizes and are packaged in high-visibility, orange, wall-mountable canisters or pouches.

BURN JEL EU

Unlike petroleum-based products, which can trap heat in the skin, Water-Jel Burn Jel EU actually draws the heat out of the burn and begins dissipating it on contact. Based on the same water-based gel technology used in our Fire Blankets and Burn Dressings, Burn Jel EU offers fast, cooling relief from the pain of minor burns. Available in unit-dose packs, squeeze bottles and a pump dispenser.



BURN KITS & BURN STATIONS

Water-Jel Burn Kits and Emergency Burn Stations are available in a variety of designs to better meet the different work environments of first responders, safety managers, and medical professionals. Developed in collaboration with paramedics around Europe, these burn kits and burn stations offer different combinations of our Burn Dressings, Burn Wrap/Fire Blankets, and Burn Jel to meet any burn emergency. All Water-Jel Burn Kits come in handy transport and storage bags. All Water-Jel Emergency Burn Stations are manufactured from durable, wipe-clean plastics.



Water-Jel Ambulance Burn Kit:

Developed in cooperation with ambulance services across Europe to match measurement proposals, the Ambulance Burn Kit supplies the rescuer with a maximized set of Water-Jel dressings and a blanket.

CONTENT: Ambulance Burn Kit		
Item	Size	Quantity
Burn Wrap	91 x 76 cm	1
Face Mask	30 x 40 cm	1
Hand Dressing	20 x 55 cm	1
Dressings	20 x 45 cm	1
Dressings	10 x 40 cm	2
Conforming Bandages		4
First-Aid Scissors		1

Water-Jel Fire Service Kit:

The right choice for first responder teams to be able to offer a first aid treatment for burn injuries. Fire rescuers are also a high-risk group for burn injuries and this kit offers a full set of burn dressings for the most common burn areas experienced with their job.

CONTENT: Fire Service Kit		
Item	Size	Quantity
Face Mask	30 x 40 cm	1
Dressings	10 x 40 cm	2
Dressings	10 x 10 cm	2
Dressings	5 x 15 cm	2
Conforming Bandages		2
First-Aid Scissors		1

Water-Jel Burn Kit XS:

The right choice to be prepared for first aid treatment of burn injuries at home. This compact, handy kit allows for easy storage in your private first aid cabinet.

CONTENT: Burn Kit XS		
Item	Size	Quanity
Dressings	5 x 15 cm	2
Dressings	10 x 10 cm	2
Burn Jel	50 mL	1
Burn Jel EU Unit Dose	3 x 4 g sachet	1



Water-Jel Emergency Burn Station - Unit 1:

Designed for use in hospital emergency rooms and departments, doctors' and health care surgeries, industrial first aid rooms, and highrisk industrial sites (e.g. sub stations).

CONTENT: Unit 1		
Item	Size	Quantity
Face Mask	30 x 40 cm	1
Dressings	10 x 40 cm	2
Dressings	10 x 10 cm	2
Burn Jel EU	50 mL pump	1
Burn Jel EU Unit Dose	4 g	5
Conforming Bandages		2
First-Aid Scissors		1

Water-Jel Emergency Burn Station - Unit 2:

Designed for use in commercial kitchens, industrial first aid rooms, and high risk industrial sites (e.g. foundries).

CONTENT: Unit 2		
Item	Size	Quantity
Burn Wrap	91 x 86 cm	1

Water-Jel Emergency Burn Station - Unit 3:

Designed for use in commercial kitchens, domestic kitchen, workshops and workstations.

CONTENT: Unit 3		
Item	Size	Quantity
Dressing	5 x 15 cm	1
Dressing	10 x 10 cm	1
Burn Jel EU Unit Dose	3 x 4 g sachet	2
Conforming Bandages		1



Unit 3 & 4, The Mead Business Centre Mead Lane Hertford, Herts SG13 7BJ United Kingdom Tel: +44 (0) 1992 583 222 Fax: +44 (0) 1992 583 229 Email: info@waterjel.net Web: www.waterjel.net

Sales Representative

Ferno	(UK) Limited		
Ferno	Ferno House		
Stubs Beck Lane			
Cleckheaton			
West Yorkshire			
BD19 4TZ			
T 1			
Tel:	01274 851 999		
Fax:	01274 851 111		

eMail: sales@ferno.co.uk Web: www.ferno.co.uk

Product Number	Size	Product Name
0206	05 x 15 cm	Burn Dressing
0404	10 x 10 cm	Burn Dressing
0416	10 x 40 cm	Burn Dressing
0818	20 x 45 cm	Burn Dressing
0820	20 x 55 cm	Hand Dressing
1216	30 x 40 cm	Face Mask
3630	91 x 76 cm	Burn Wrap
P3630	91 x 76 cm	Burn Wrap (Pouch)
7260	183 x 152 cm	Fire Blanket-Plus
P7260	183 x 152 cm	Fire Blanket-Plus (Pouch)
9672	244 x 183 cm	Heat Shield
BJ50	50 mL	Burn Jel Squeeze Bottle
BJ120	120 mL	Burn Jel Squeeze Bottle
BJ240	240 mL	Burn Jel Pump Bottle
BJ340	3 x 4 g	Burn Jel Unit Dose
BJ640	6 x 4 g	Burn Jel Unit Dose
BK-LG-WJE		Ambulance Burn Kit – In Pouch Content: 1 x P3630, 1 x 0818, 1 x 0820, 2 x 0416, 1 x 1216, 4 x SCB1, 1 x TCS1
BK-SM-WJE		Fire Service Burn Kit – In Pouch Content: 2 x 0206, 2 x 0404, 2 x 0416, 1 x 1216, 2 x SCB1, 1 x TCS1
BK-XS-WJE		Burn Kit XS – In Pouch Content: 2 x 0206, 2 x 0404, 1 x BJ50, 1 x BJ340
BK-SM-IBK		Industry Burn Kit – In Pouch Content: 2 x 0206, 2 x 0404, 1 x 0416, 1 x 1216, 1 x BJ50, 1 x SCB1, 1 x TCS1
BK-SM-CBK		Catering Burn Kit – In Pouch Content: 2 x 0206, 2 x 0404, 1 x 0416, 1 x BJ50, 1 x SCB1, 1 x TCS1
BK-SM-TBK		Water Jel Training Burn Kit – In Pouch (Nonsterile, for education only)
BKCK-SM-WJE		Critical Burn Kit – In Red Pouch Content: 2 x 0820, 2 x 0404, 1 x 1216, 2 x SCB1, 1 x TCS1
EBS-Unit1		Emergency Burn Station – In Wall Station Content: 2 x 0404, 2 x 0416, 1 x 1216, 1 x BJ50, 5 x BJ140, 2 x SCB1, 1 x SS1
EBS-Unit2		Blanket Dispenser Station – In Wall Station Content: 1 x P3630
EBS-Unit3		Work Station – In Wall Station Content: 1 x 0206,1 x 0404, 2 x BJ340, 1 x SCB2

Rev# wj-cat-002-2