



EZ-GLIDE LBS POWERTRAXX Bariatric Powered Stair Chair



The **NEW** industry-leading patient handling of the Powered EZ-Glide LBS tracked chair... now in a stronger, wider specification for larger patients

Moving patients safely up and down stairs can be distressing to the patient and carry substantial manual handling injury risks. Ferno has the product to aid in these problems, also offering more dignity and comfort when moving the patient – the EZ-Glide LBS PowerTraxx™.

EZ-Glide LBS PowerTraxx™ is so much more flexible than any other powered tracked chair. Capable of moving down and upstairs under controlled power, you can slow it down... turn it around... bring it to a halt... all at the touch of a button.



- **Effortless patient handling** – adjustable speed control for safe, easy use; touch-button wireless control panel for convenient operator use
- **Stronger and wider** – incorporating extended sides, a reinforced foot-plate and a strengthened fixed frame for patients with a greater physical mass
- **Uniquely easy-handling** – 25 mm track-to-ground clearance and smooth 100 mm front swivel wheels
- **High load capacity** – will transport patients and equipment loads of up to 227 kgs (35 stone)
- **Convenient to use** – detachable battery stays charged over 20 full flights of steps, just 60 minutes recharge time
- **Excellent patient safety** – an automatic braking system combined with options such as adjustable cushioned headrest/shoulder support and extended footrest, ensures patient safety and comfort
- **Flexible and versatile** – climbs and descends stairs, goes around corners, tackles difficult terrain... entirely without strain

EZ-GLIDE LBS POWERTRAXX

Bariatric Powered Stair Chair



Features & Benefits

- Easily manoeuvrable in tight spaces – 150 mm rear wheels for effortless handling even on most stair landings and other confined areas
- Reassuring patient security – with Biosafe 4-point harness and extending footplate complete with straps
- Feature rich – including rear folding handles, extended footplate, IV pole and winch attachment, all included as standard
- Complete operator control – easy to operate powered track provide ease of handling
- Improved comfort – supplied with seat and backrest cushion, and padded headpad. Headrest/shoulder support also available
- Multi-purpose – fixed frame chair, can be used as a in-vehicle patient seat when used in conjunction with NMI Saferider* restraint system and Pocket Lock* or Floor Anchor Strap System*

Specifications

HEIGHT (max)	950 mm
HEIGHT (min)	610 mm
DEPTH	850 mm
WIDTH	595 mm
WEIGHT	37.9 kg
PATIENT WEIGHT CAPACITY	227 kg / 35.7st / 500 lbs
BATTERY TYPE	28V Lithium Ion

Accessories



- Milwaukee battery 5Ah
75-0710-001



- 240v mains Milwaukee battery charger
75-0711-002



- Saferider restraint & Pocket lock
SAF-B-P



- Chair floor bracket c/w Pocket Lock
PL-FEGCS

In-Vehicle Fixing Options

Part Numbers

073132940	EZ-Glide LBS with PowerTraxx, Red fixed-frame, non-folding, c/w 4-point Biosafe restraint straps and extended footplate
031402600	Replacement cushion for EZ-Glide LBS PowerTraxx
031385825	Replacement headrest for EZ-Glide LBS PowerTraxx
038300000	Dust cover for EZ-Glide LBS PowerTraxx
082197700	D-size cylinder Oxy-Clip2 Bottle Holder
082197900	CD-size cylinder Oxy-Clip2 Bottle Holder
65-2004-002	Replacement connect cable for EZ-Glide PowerTraxx
65-1406-001	Replacement pouch to hold connect cable
65-1400-024	Replacement wireless antenna module for PowerTraxx