



**Wheel-Kit 71 W for
FERNO Basket Stretcher 71/71 S**

WHEEL-KIT 71 W



**Wheel-Kit 71 W for
FERNO Basket Stretcher 71/71 S**

The wheel kit for our basket stretcher 71 or 71 S facilitates transport over longer distances.

With a few simple steps, the kit is attached to the foot part of the basket stretcher with 4 quick-release fasteners.

The 260 mm big wheels with all-terrain profile and pneumatic tires ensure easy running even on uneven terrain.

Powder coated aluminum frame.



GENERAL SPECIFICATIONS

	Length	Width	Height	Weight	Load Capacity
Wheel-Kit 71 W	60 mm	60 mm	32 mm	4.5 kg	200 kg