



Fiber Optic Laryngoscopes

FERNO2 / Spiracle Technology laryngoscope blades and handles are designed with medical professionals in mind. Constructed of stainless steel, these durable fiber optic laryngoscope blades are of the highest quality. All of the FERNO2 / Spiracle Technology laryngoscope blades are compatible with ISO “Green System” handles.

- FERNO2 / Spiracle Technology laryngoscope blades come in a wide variety of Mac or Miller sizes
- Compatible with ISO “Green System” Handles
- The bulb-in-handle design provides the user with a long-lasting, bright light source while removing the potential of bulb aspiration found in “bulb on blade” configuration models.

Baseline / General Specifications

Ergonomic shape/design. High transparency bright fiber optic bundle, “no” bulb at distal end of blade. XENON Halogen bulb in handle. Helps reduce bulb aspiration and potential hot spots/burns. Compatible with ISO 7376 (Green System) handles.

Material: Stainless Steel (satin finish for low glare), Fiber Optic bundles- fused silica.

Available Blade Sizes Part Number

MacIntosh 0	T22427-0	
MacIntosh 1	T22427-1	
MacIntosh 2	T22427-2	
MacIntosh 3	T22427-3	Overall Length 126mm, Inside Length 101mm
MacIntosh 4	T22427-4	Overall Length 160mm, Inside Length 135mm
Miller 0	T22428-0	
Miller 1	T22428-1	
Miller 2	T22428-2	
Miller 3	T22428-3	



A FERNO Group Company

10601 Calle Lee, #190

Los Alamitos, CA 90720

Phone: (714) 418-1091 • Fax: (714) 418-1095

e-mail: info@spiracle.com • www.spiracle.com