

# FERNO SPED II

## Spinal Protection Extrication Device

DEVELOPED and DESIGNED for **EXTRICATION** from vehicles and other confined spaces



FOR MORE DETAILED  
INFORMATION SEE OVER

### Specifications

**FERNO** – is a world leading innovator and manufacturer of patient transfer and rescue solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **SPED** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

	Length (packed)	Width folded	Width unfolded	Weight	Load Limit
SPED II	1000 mm	650 mm	670 mm	2.7 kg	200 kg/31st/440lbs
SPED II in Bag	700 mm	250mm (depth)	300 mm	4.7 kg	

**call** +44 (0) 1274 851 999  
**email** sales@ferno.co.uk  
**www** ferno.co.uk

 **FERNO**<sup>®</sup>

# FERNO SPED II

## Spinal Protection Extrication Device

For over 40 years the KED set the standard for the extrication device to use for suspected vehicle spinal and neck injury accidents. From this original concept we have created SPED, an amazing leap forward in extrication safety, incorporating vacuum mattress immobilisation technology.

- **Low Friction Material** - easier deployment when introduced in between the patient and seat
- **Colour Coded Strapping** - buckles and straps stored within the device ensuring easy deployment by avoiding tangling of loose straps
- **Vacuum Immobilisation Technology** - within seconds the device becomes rigid, moulding to patients' contours
- **Integral Spinal Stiffening Ribs** - gives maximum support during extrication
- **Integrated Head Immobilisation** - superior C-spine immobilisation helps reduce risk of further injury whilst moving the patient
- **Seven Lifting Handles** - provides more points of access for better lifting and manoeuvrability



### Part Number

23005001001

Ferno SPED II, "ReadytoGo" complete with introducer, restraint strap re-locator, repair kit + carry bag

FERNO - SPED II FPC002 10.15

**call** +44 (0) 1274 851 999  
**email** sales@ferno.co.uk  
**www** ferno.co.uk



Ferno (UK) Limited, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.