

Ferno CCT-PX Critical Care Transport Trolley



Specialist Critical Care Transfer Trolley made to your specifications

CCT-PX is Ferno's latest enhanced Critical Care trolley, designed and built to individual requirements but configurable to ever-changing clinical needs, ensuring that all the medical equipment required for each patient transfer is safely attached to the trolley at all times. Feature-rich in design, your bespoke custom-made CCT-PX can be transformed in real-time according to the immediate needs of each transfer.

The trolley can even be configured to incorporate Bariatric capabilities, using the optional SX Surface eXtender cotsides.

Specialist equipment such as Balloon Pumps, ECMO machines and cooling/heating systems can also be accommodated at the time of design.



PacRac+ equipment table



Side-Rail Interface



Ferno CCT-PX Critical Care Transport Trolley

Overview

- Designed to Specification Bespokemade according to each individual's requirements, but with the flexibility to accommodate varying needs
- Integrated Infection Control clean lines, wipe-clean pressure-reducing ABPR mattress, biosafe patient restraints and anti-bacterial paint surfaces all contribute to a high level of infection control and safety
- Quick Configurability integrated trolley side-rail interface allows a host of options and accessories to be easily attached to the trolley
- CEN Tested Safety dynamically tested and certified to EN 1789, compatible with Ferno's complete range of 2-Part Locking Systems

Specifications

LENGTH (fixed)		1970 mm
HEIGHT (fixed)	To the underside of the mattress	825 mm
WIDTH	Normal	570 mm
	SX Cotsides for Bariatrics	920 mm
Load Limit (SWL)	Patient*	250 kg / 40 st
WEIGHT	Mattress and restraint not included	69 kg**

Ferno's CCT-PX comes with the most advanced range of features and configurable options ever developed in a Critical Care Patient Transfer trolley

CCT-PX Baseline Features

- · Pole mounts to hold medical equipment
- Side-interface for attachment of accessories
- Olson 12-gang Power Distribution Unit
- 2 * E-size cylinder holders (with CD adapter)
- Lateral-folding cotsides (option of SX cotsides)
- Black webbing 'Cross-strap' patient harness



Vertical Fixing Poles For attachment of pole-clamped medical devices



Adjustable Winged Headrest

Providing comfort and adjustable to different size patients



E-size Cylinder Holders Holds E size cylinders and also CD size with supplied

Anti-Bacterial Pressure-Reducing mattress

Detachable adjustable winged headrest

Adjustable backrest & leg position platform

(Trendelenburg, Shock & Knee Contour)

• Stowable foot-end Push-poles (2)

Brakes on all 4 wheels



adapter kit



Power Distribution Unit Easy and safe supply of power for all the medical

Anti-static wheel

· Bump wheels on all 4 corners

*Maximum equipment weight dependant on trolley configuration * Standard mass (weight dependant on final configured specification)

- Backrest monitor hook
- · Winch connector for loading in vehicle



Lateral-Folding Cotsides Simply fold down when transferring a patient



Trolley side-interface

Allows various options and accessories to be attached to the trolley, as and when required

Configurable Options & Accessories

- Wipe-clean Biosafe Cross-Strap harness .
- SX Surface Xtender ratcheting cotsides with extendable patient surface
- Clip-on Spilt-Scoop Extreme Bariatric surface
- Detachable 3-stage telescopic IV Pole



SX Surface Extender Provides added comfort and security for patients of any size FWESX



3-Stage IV Pole Quickly place an IV bag anywhere around the patient LAIV3SCL

- Dedicated O2/ZX/E/Nitric cylinder holder
- PacRac+ equipment table
- Wide range of medical device brackets, fixed & swivel mount
- · Drawer unit / cabinet options





CD Cylinder Holder Dedicated CD or Nitric Oxide cylinder holder attached to trolley baseplate

- Head-end stowage net
- Additional 2 Push-Poles for head-end
- Side-fixing manoeuvring handles
- · Directional wheel locks at head-end
- Balloon Pump & ECMO options



PacRac+ Equipment Table Attach monitoring devices & other equipment, opens both sides FWEPR+IT



Swivel Mounts for PacRac+ Allows equipment to be mounted with 360 degree rotation allowing devices to be viewed from various positions

