



POWERX™

 **FERNO®**

POWERX™ EXTRA POWER



Designed for

**ENHANCED
MANOEUVRABILITY**

With improved stretcher weight distribution and wheel configuration, in combination with features such as ergonomic push-poles and directional wheel locks, making **POWERX** easier to manoeuvre

Configurable to give

**FULL BODY BARIATRIC
CAPABILITY**

Offering ratcheting Surface eXtending cotsides and Split-Scoop patient surface extension options, combined with a 350kg (55 stone) unassisted-lift trolley SWL, helping to reduce back injuries, safeguarding carers, retaining staff and saving costs

Improved stretcher weight distribution and

**FLAT PATIENT
SURFACE**

Non-contoured flat Anti-Bacterial Pressure Reducing mattress provides the optimum stretcher height, making lateral patient transfers fast, safe and easy



Anti-bump wheels on all four corners



Detachable 28V DC, 5Ah Battery

**ONE SOLUTION FOR
ALL PATIENT JOURNEYS**

FOR MORE INFORMATION



FERNO

POWERX

FERNO



Wheel brakes on the frame, making them easy to access. Directional wheel locks on both head-end wheels

INTERCHANGEABILITY ACROSS YOUR FLEET

Designed and built with

LESS COMPONENTS

Increasing reliability, combined with 12 month service intervals and minimum replacement part requirements reduces equipment whole-life cost



Latching defibr hook for safe handling of equipment



Side-lighting for improved visibility



Easy to operate backrest positioning



Stowable Push-Poles

BATTERY POWERED TECHNOLOGY

5Ah Lithium-Ion battery - quick to recharge and long-lasting both during daily operation and throughout the stretchers life



FOR MORE INFORMATION



ferno.co.uk | +44 (0) 1274 854 510 | sales@ferno.co.uk

POWERX IS THE MOST ADVANCED POWERED X-FRAME STRETCHER ON THE MARKET, DEVELOPED TO ENHANCE PATIENT CARE



Storable head and foot-end Push-Poles



Directional wheel locks at head-end



Manual release handle



Winch-link for assisted loading on ramps

POWERX is designed and built in the UK



***Adjustable winged headrest**



***Backrest monitor hook**



***CD/ZA O2 cylinder holder**



***Storage net head-end**



***Split-Scoop Large Patient Surface**



***Detachable 3-stage telescopic IV Pole**



***PacRac+, for holding medical equipment**



***Side-fixing manoeuvring handles**



***12v or 230/240v battery charger options**




***Lateral-folding cotside**



***SX Surface Xtenders ratcheting cotside (Extended patient surface)**



***Side-stretcher lighting**



***Electrically-operated backrest and trendelenburg platforms**

ACCESSORY SYSTEM

PacRac™ Equipment Mount
Medical devices are safely secured - the opening side-supports allow for patient transfers without removing the PacRac+

SX Surface eXtenders
Provides added comfort and security, accommodating patients of varying shapes and sizes

3-Stage IV Pole
Allows you to quickly place an IV line and bag anywhere around the patient using the universal accessory rail

Adjustable Headrest
Designed for improved lateral transfers and enhanced patient comfort

Latching Equipment Hook
Hang your monitor on the backrest and let **POWERX** carry the weight for you

Storage Net
Provides additional storage for added convenience

Integrated O² Holder
Securely attaches an oxygen cylinder to **POWERX**

Universal Accessory Rail
Allows you to move and attach accessories anywhere along the sides of **POWERX**



Specifications

LENGTH	Short to Full	1450-1930 mm
HEIGHT	Low to Raised	480-880 mm
WIDTH	Normal	570 mm
	SX Cotsides Horizontal	920 mm
	Split Scoop kit attached	945 mm
Load Limit (SWL)	Raised Limit	350 kg / 55 st
	Lowered Limit	400 kg / 63 st
BATTERY	28 Volt DC	5Ah
WEIGHT	148 lbs	78 kg*
FASTENING SYSTEM	Two Part Lock	13 kg*

* Weight dependant on final configured specification. POWERX illustrated shown with various optional configurations and accessories. For full available range see configurable options and accessories list.

POWERX has been designed to meet European dynamic test standards - increasing safety for you and your patient. CEN EN 1789 certified

Fully supported Ferno Care service and maintenance with a full range of tailor-made packages available

DISCLAIMER. * Optional accessories available. Pictures of products may differ. Contact Ferno UK for the latest features and available accessories. Ferno reserves the right to change product specifications without notice.

FOR MORE INFORMATION



ferno.co.uk | +44 (0) 1274 854 510 | sales@ferno.co.uk

UK 05/19 - R2