FERNO°



The **POWER X1 Ambulance Cot** offers more power, stability, maneuverability, and speed than the leading competitive powered cot.

With an unassisted lift capacity of 700 lb, the POWER X1 is the ultimate solution for both everyday and bariatric transport. Its industry-leading battery technology provides high-powered, reliable performance on every call.

The POWER X1 is cross-compatible, as it can be used with a range of fastener systems.



ENDURACharge Power System delivers 113 cycles* per charge. *One cycle equals one full lift and one full lower



RAPIDLift Retraction System achieves one full cycle in just 3.4 seconds.



ERGOFrame Flexible Positioning shortens the frame by as much as 22 in.



82% more STABILITY



ADVANCED STABILITY TECHNOLOGY

40% more SPEED



9% more MANEUVERABILITY



**Stryker Power-PRO™XT











Preventative maintenance program - Parts, labor, travel, two annual preventative maintenance inspections, unscheduled service, and product health report. Replacement parts do not include: mattresses, batteries,** accesories, other disposable or expendable parts, or damage due to misuse. **Batteries are included if accompanying Battery Plan is purchased.



Ambulance Cot





FEATURES AND BENEFITS

ENHANCED PATIENT CARE, COMFORT AND SAFETY

- Pneumatic backrest adjusts to a comfortable position for each patient
- Infinite cot-height adjustment gives you greater flexibility
- Shock frame allows you to raise patient's feet for treatment of shock
- Knee gatch (option) increases patient comfort and clearance between foot-end operator and patient's feet
- Sealed bolster mattress is designed for patient comfort and features a sanitary Hook and Loopfree design
- SAE-compliant FERNO X-Restraints are color-coded for easy application

CROSS-COMPATIBLE DESIGN

PRO F1 Compatibility Kit allows you
to secure the POWER X1 to the PRO
F1 Universal Cot Fastener and Stryker
Performance-LOAD™. The POWER X1
is also compatible with FERNO and
Stryker antler-and-rail fasteners.

IMPROVED HANDLING AND MANEUVERABILITY

• 5-position articulating ERGOFrame enables optimal ergonomic positioning for operators

- 6-by-2-inch oversized wheels with wheel locks on the control end and optional directional wheel locks on the head end offer maneuverability, even over rough terrain
- Control-end wheel locks prevent cot rolling during patient transfer and medical procedures
- Two sets of contoured, ergonomic controls and Sure Grip™ dipped lift bars and grasping surfaces enable easy handling
- Backrest-mounted equipment hook helps keep your hands free

BUILT FOR LONGEVITY AND SIMPLE OPERATION

- Tested in a controlled environment to 10+ years of service
- Durable frame is made of powdercoated aluminum
- Programmable load height adjusts to match the height of your ambulance floor
- ICS® Integrated Charging System (available with Pro F1 Universal Cot Fastener) eliminates battery management
- Innovative Cycle Tracker™ records usage to help you keep up with maintenance and extend equipment life
- Manual back up mode provides complete operation if power loss occurs

SPECIFICATIONS

Item Number	Base model: 0015807
Length (max)	80 in / 2030 mm
Length (min)	58 in / 1483 mm
Weight*	143 lb / 65 kg
Weight Capacity	700 lb / 318 kg (Unassisted) 1,100 lb / 498 kg (In collapsed position)
Power	36 V Lithium-lon

* Cot is weighed without mattress and restraints.

STANDARD FEATURES

- Equipment Hook
- Bolster Mattress
- Shock Frame
- Foot-end Wheel Locks
- Foot-end Pull Handle
- FERNO X-Restraints
- Fold-down Sidearms
- · Antler and Rail Compatibility

COMPLIANCE

The **POWER X1** has been designed to meet international test standards to increase safety for you and your patient.

- SAE J3027
- CAAS-GVS v1.0
- IP56
- KKK-A-1822F

- NFPA 1917
- AMD 004
- IEC 60601-1-2
- ISTA 1B

The POWER~X1~is~compliant~with~the~above~standards~when~used~with~the~PRO~F1~Compatibility~Kit.

OPTIONAL ACCESSORIES



Three-Stage IV Pole No. 1590160

SX Surface Xtender™

Sto-Net[™] Load Frame

No. 1590159

No. 0822482



PacRac*+ Equip Mount
No. 1590161



Replacement Battery No. 0822484



Side Lift Handles No. kit 0821998



Head-End Oxygen Bottle Holder No. 0822485



Cycle Tracker No. 0832159



Directional Wheel Locks No. 0900632

BARIATRIC ACCESSORIES



LBS Jr. No. 0822499



Mattress LBS Jr. No. 0314050



Model 513- LBS IV Pole No. 0087173

ACCESSORIES NOT PICTURED: PRO F1 Compatibility Kit, Patient Shield & Knee Gatch

Power-PRO XT and Performance-LOAD are trademarked products of Stryker Corporation FERNO is not associated with Stryker Corporation.