

FERNO Paraguard Shoulder Harness

For use with Paraguard stretcher handles to give better support
& reduce lifting effort required



Paraguard shoulder harness is designed for use with Paraguard handles to aid rescuers, by distributing the weight being carried, reducing the lifting effort required and enhancing the rescuers' body positioning. The harness is simply worn around the upper body, clips located at the ends of the incorporated lifting straps are attached to dedicated fixing points on the ends of the Paraguard carry handles. The lift straps length can be adjusted according to required positioning and comfort.

010772801	Paraguard Shoulder Harness
Weight	1.3 kg (2 harnesses)
Load rating	136 kgs / 21 stone / 300 lbs
Dimensions	56 cm * 32 cm (body panel dims)

FERNO - Paraguard Shoulder Harness FPC001 01.21

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

