

FERNO Neomate[®]

Paediatric Harness System

Transport infants safely and securely on a trolley in our easy-to-use harness system



FOR MORE DETAILED
INFORMATION SEE OVER

Specifications

Length	Width	Height	Weight	Load Limit
533 mm (open) 330 mm (folded)	266 mm (open)	88 mm	1.1 kg	2.3 kg (7 lb) Min 6.4 kg (14 lb) Max

FERNO – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **Neomate[®] Harness** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER.**

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



FERNO Neomate® Paediatric Harness System

Safe, comfortable and versatile

Fits securely onto a wide range of Ferno trolleys

Safety is paramount during the emergency transportation of infants. Your chosen harness needs to support infants securely and offer the versatility to accommodate a broad range of weight and sizes.

The **Neomate®** harness is lightweight and fully adjustable, making it easy to use. It is also comfortable for the patient, ensuring that the infant is not upset by feeling overly restricted.

Neomate® is designed to accommodate infants weighing between 2.3 kg and 6.4 kg, (or 7 to 14 lbs). It is easy to use and made from durable, easy clean vinyl and folds in half for easy storage.

As you would expect it is compatible with ALL patient transfer trolleys.

- **Simple and Secure** –
five-point fully adjustable harness system clips snugly in place
- **Quick and easy to use** –
central easy-release button
- **Lightweight** –
adjustable padding is made from soft, lightweight vinyl
- **Safe and hygienic** –
made from non-toxic, easy-clean materials
- **Tested on Ferno trolleys** –
to meet CEN dynamic crash test standards



Side padding for smaller baby

Part Nos:

031-4114CF

Neomate® Paediatric Harness System

MANU000140

iNX adapter bracket for using Pedi-Mate, Pedi-Mate PLUS & Neomate on iNX



FERNO - Neomate® FPC006 08.21

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO®**

Ferno (UK) Limited, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ, England