

FERNO EasyFix II Vacuum Mattress

Revolutionary body shape sets new standards in
Full-body Vacuum Mattress immobilisation



**FOR MORE DETAILED
INFORMATION SEE OVER**

Specifications

Length	Width	Weight	Load Limit
2080 mm	1300 mm (Max)	8.5 kg (incl. straps) 9.6 kg (incl. straps, footrest & pump)	250 kg/39st/550lbs

FERNO – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **EasyFix Vacuum Mattress** typifies our passion for top quality, highly functional design...

patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO**

FERNO EasyFix II Vacuum Mattress

Achieve secure immobilisation with improved vacuum mattress design

EasyFix ii utilises advanced technological developments to deliver superior immobilisation, exactly where and when it's needed. The system uses a V-shaped body design to optimise the vacuum mattress capabilities, whilst at the same time drawing on enhanced internal multi-chamber technology for superior stability throughout. The revolutionary X-shape restraint system immobilises the upper body, hips, pelvis and legs, demonstrating the versatility of this unique mattress. An enhanced angled headstrap on the EasyFix Mk II model makes deployment even easier and more effective.

However sophisticated an immobilisation technique may be, it also needs to be quick and easy to deploy. **EasyFix** offers improved performance here too, as the suction pump device fits directly onto the mattress valve so no time is wasted attaching adapters – and easily adjusted strap system ensure operators can save time when it's needed most.

- **Superior immobilisation** – unique chamber design and body shape with built-in wooden upper support to hold the patient in the correct position more securely. Shaped head and shoulder support also provides further patient safety
- **Improved immobilisation of upper back, neck and head** – reinforced C-spine base plate
- **Durable and highly effective** – double layered material helps maintain vacuum longer and reduces accidental damage or punctures
- **High carrying capacity** – load limit 250kg
- **X-ray transparent, MRI compatible** – no need to move the patient out of the vacuum mattress for scanning
- **Easy storage** – folds easily to take up 20% less space than other designs of vacuum mattress
- **Simple maintenance** – Polyester fibre, PVC coated
- **Smooth handling** – with 10 ergonomically designed handles



Part Nos:	
271501100	EasyFix II Vacuum Mattress
271501101	EasyFix ii Vacuum Mattress Kit (includes footrest pump, repair and bag)

FERNO - EasyFix II Vacuum Mattress FPC001 03.21

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



Ferno (UK) Limited, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.