



## Integrated Patient Transport & Loading System



- **Dual, powered actuators** maximize tip stability and lift patient loads up to 700 lb
- **Simple one-touch operation** helps you stay focused on patient care
- **SAE J3027 compliant** when you combine the iN/x with the iN/JLINE Fastening System
- **Telescoping load frame** maximizes maneuverability in tight spaces
- **Integrated Charging System** eliminates the need for battery management

Everyday operations in EMS — like transferring patients or loading and unloading cots from an ambulance — put tremendous strain on the backs of EMS professionals.

The **iN/x Integrated Patient Transport & Loading System** is designed to enhance the safety of both patients and EMS professionals by effectively eliminating lifting. EMS professionals can lift, lower, load, and unload with the push of a button, taking the strain off their backs.



### No-Lift Design

allows you to load and unload with ease



### Independent Leg X-Frame Design

allows you to navigate obstacles with your thumb and not your back



### Integrated Lighting System

with continuous and flashing modes allows you to see and be seen in dark and hazardous conditions.

Certified  
SAE J3027  
Compliant



Preventative  
Maintenance  
Program



Preventative maintenance program - Parts, labor, travel, two annual preventative maintenance inspections, unscheduled service, and product health report. Replacement parts do not include: mattresses, batteries,\*\* accessories, other disposable or expendable parts, or damage due to misuse.

\*\*Batteries are included if accompanying Battery Plan is purchased.



## Integrated Patient Transport & Loading System



### FEATURES AND BENEFITS

#### POWERFUL, STABLE PATIENT HANDLING

- **Dual powered actuators** maximize tip-stability and can lift patient loads up to 700 lb.
- **Automatic high-speed extend & retract** reduces load and unload time
- **Chair position** facilitates transfer of seated patients and COPD patients
- **Innovative X-Frame design** simplifies positioning with just the push of a button

#### PATIENT COMFORT AND SECURITY

- **Patient surface** is level at all heights to enhance procedural care
- **Pneumatic-assisted backrest** provides patient comfort
- **Shock position** raises patient's feet for treatment of shock

#### TIGHT SPACE MANEUVERABILITY

- **Telescoping load frame** maximizes maneuverability in tight spaces

- **360° locking swivel wheels** allow movement in any direction at any bed height
- **Telescoping foot-end handles** increase patient foot clearance

#### SIMPLE, SAFE, SECURE

- **Simple one-touch operation** helps you stay focused on patient care
- **Programmable load height** adjusts to match your ambulance floor height
- **ICS® Integrated Charging System** eliminates need for battery management
- **LCD Display** with battery charge indicator provides system status and operational cues
- **Continuous and Flashing Lights** allow you to see and be seen in dark and hazardous conditions.
- **In-fastener shut-off** prevents accidental operation in the ambulance
- **Memory-free battery and manual back-up system** keep the iN/X ready for every call

### SPECIFICATIONS

<b>Item Number</b>	Base model: 0000INXINLINDWL
<b>Length (max)</b>	81 in / 2057 mm
<b>Length (min)</b>	63 in / 1600 mm
<b>Weight*</b>	203 lb / 92 kg
<b>Weight Capacity</b>	700 lb / 318 kg (Unassisted)
<b>Backrest Adjustment</b>	0° — 76°
<b>Shock Frame Adjustment</b>	0° — 16°
<b>Power</b>	36 V

\* Cot is weighed without mattress and restraints.

### COMPLIANCE

The iN/X has been designed to meet international test standards to increase safety for you and your patient.

- BS EN-1789
- AS/NZS-4535
- IP56
- SAE J3027
- KKK-A-1822
- NFPA-1917
- IEC 60601-1
- CAAS GVS v1.0

### OPTIONAL ACCESSORIES



**Three-Stage IV Pole**  
No. 1590160



**PacRac+ Equip Mount**  
No. FWEPRIT



**Directional Wheel Lock**  
No. FWECDL



**SX Surface Xtender™**  
No. 1590159



**Equipment Mount**  
No. 0822451



**Load Frame Handles**  
Side No. FWESLH  
Front No. FWELFH



**Sto-Net™**  
No. 0822453



**Patient Shield**  
No. 0822097

### BARIATRIC ACCESSORIES



**LBS Jr.**  
No. 000LBSJRINX



**Mattress LBS Jr.**  
No. 0314050



**Model 513- LBS IV Pole**  
No. 0087173

ACCESSORIES NOT PICTURED: Load-End & Operator End Leg Storage Nets