TRANSCEND™ + POWERTRAXX®



Empowering every step.

The Transcend™ Stair Chair with PowerTraxx® provides maximum safety and power when transporting a patient up or down stairs.
It's electronically-controlled motor, direct drive transmission, and powered tracks push and pull patient loads of up to 550 pounds. The chair and track system carry the weight, increasing patient safety and reducing oeprator fatigue and risk of injury.

POWERCLIMB

TRANSPORT AN AVERAGE OF 40 FLIGHTS OF STAIRS WITH A TOUCH OF A BUTTON

DESCENTCONTROL

LIMITS TRACK SPEED DURING DESCENT OR ASCENT TO A SAFE, MANAGEABLE LEVEL

EZGLIDE®

PROVIDED COMFORTABLE RIDE ON ALMOST ANY SURFACE WITH 60% LESS SHOCK THAN COMPETITIVE STAIR CHAIR

STAIRBRIGHT

INTEGRATED LED SCENE LIGHTS TO ILLUMINATE SURROUNDINGS UP TO 1 FLIGHT OF STAIRS WITH TWO BRIGHTNESS OPTIONS.

FERNO

TRANSCEND™ + POWERTRAXX®

ENHANCED PATIENT CARE, COMFORT, AND SAFETY

- Extending footrest and lower track angle increase patient scense of security
- Patient seating surface is contoured and has removable injection molded seat panels
- Patient restraints for security and safety
- Armrests offer additional support and security for patients while enchancing comfort and positioning during transport or transfers*
- Headrest provides added comfort and stability fo rthe patient's head and neck during transport, helping reduce movement and increase confidence*

INCREASED PARAMEDIC SAFETY, EFFICIENCY, AND MANEUVERABILITY

- Smooth 5-inch front swivel wheels and extra tall 8-inch rear wheels with integrated shock absorbing tires that navigate rough terrain and do not collect debris.
- 5-position lift bar and handles reduce the strain of lifting and provide maximum control for operators of various heights (wireless motor controls work in any handle position)*
- built-in electronic speed control limits track speed during descent to a safe, manageable level
- One-inch track-to-ground clearance provides smooth rolling over carpet and rough surfaces, as well as easier descent down stairs
- Easy to operate with push-button wireless controller
- Rear lift handles provide improved ergonomics for safer, more comfortable lifting and greater control*
- Posititve locking hinge mechanism ensures chair is locked in place whether deployed for use or folded for carrying and storage

ENGINEERED FOR RELIABILITY, DURABILITY, AND PERFORMANCE

- Two battery options: ENDURACharge 36v runs up to 40 flights of stiars per charge or the Milwaukee® 28v lithium-ion battery that is commercially available and runs up to 30 flights per charge.
- All-metal frame constructikon provides durability and long life
- 2-year warranty

*OPTIONAL ACCESSORIES

SPECIFICATIONS

ITEM NUMBER*
0731405

HEIGHT (MAX)
60" / 1518 MM

HEIGHT (MIN)
36" / 908 MM

WIDTH (OVERALL)
20" / 505 MM

WIDTH (SEAT)
16" / 408 MM

DEPTH (FRONT TO BACK)
30" / 760 MM

LENGTH (FOLDED)
36" / 912 MM

WIDTH (FOLDED)
20" / 507 MM

DEPTH (FOLDED)
7.9" / 200 MM

WEIGHT
51 LB / 23.3 KG

LOAD CAPACITY
550 LB / 250 KG

DEPTH (FRONT TO BACK)
30" / 760 MM

POWER CONFIGURATION OPTIONS
28-VOLT DC LITHIUM-ION MILWAUKEE BATTERY
36.3 VOLT LiFePO4 ENDURACHARGE BATTERY

*STANDARD BUILD AND EQUIPMENT SPECIFICATIONS APPLY TO ITEM NO. 0731405.

ADDITIONAL BUILD OPTIONS AVAILABLE UPON REQUEST.

FERNO® RESERVES THE RIGHT TO CHANGE SPECIFICATIONS WITHOUT NOTICE.

OPTIONAL ACCESSORIES

ARMRESTS
HEADREST
O2 HOLDER
ENDURACHARGE BATTERY
ENDURACHARGE CHARGER
M28 BATTERY
M28 CHARGER

STANDARD EQUIPMENT*

REAR WHEEL LOCKS
FLIP-UP HANDLES
TELESCOPIC U-BAR WITH R&L IV CLIP
FOOTREST
GUARD BAR
M28 US BATTERY & CHARGER
RESTRAINTS

© 2025 FERNO | WWW.FERNO.COM | U.S. 1-877-733-0911 | INTL. 1-937-382-1451 | OCT-2825-V8