BARIATRIC
SOLUTIONS




Safe transport shouldn’t have limits.

Large patients require added support for their larger mass to provide greater comfort and safe, stable transfers and transport with your cot.

FERNO Bariatric Solutions allow you to quickly and easily convert standard FERNO cots into bariatric transport solution-on-scene, eliminating the need to buy a separate, dedicated stretcher.

FERNO

BARIATRIC SOLUTIONS

	iNX	POWER X1	POWER X2	35X PROFlexx
				
	No. 000LBSJRINX | No. 0015819 (with SX) No. 009008250 (with Scoop Stretcher)	No. 000LBSJRPX1 No. 0015819 (with SX)	No. 000LBSJRPX2 No. 0015819 (with SX)	No. 0822184
Standard Width	24 in / 609 mm	24 in / 609 mm	24 in / 609 mm	24 in / 609 mm
Solution	 or  or SX Surface Xtender™ LBS Jr.™ OR 65XL Split-Scoop Stretcher	 or SX Surface Xtender™ OR LBS Jr.™	 or SX Surface Xtender™ OR LBS Jr.™	 LBS Jr.™
	 			
Width (Open)	35 in / 889 mm (SX) 38 in / 965 mm (Scoop) 41 in / 1040 mm (LBS JR)	37 in / 939 mm (SX) 41 in / 1040 mm (LBS JR)	37 in / 939 mm (SX) 41 in / 1040 mm (LBS JR)	41 in / 1040 mm (LBS JR)
Width (Folded)	24 in / 609 mm (SX) 29 in / 749 mm (LBS JR)	24 in / 609 mm (SX) 29 in / 749 mm (LBS JR)	24 in / 609 mm (SX) 29 in / 749 mm (LBS JR)	29 in / 740 mm (LBS JR)

FEATURES & BENEFITS

The SX Surface Xtender™ offers 5 locked positions, and widens the cot surface of the iNX, POWER X1, and POWER X2.

The LBS Jr.™ provides a narrow profile and folding side frames for greater maneuverability through doorways.

The **LBS JR Bariatric Board** fold for easy storage in an ambulance when not in use.**iNX, POWER X2, POWER X1, and 35X PROFlexx.**

Features 1,100 lb load capacity in the lowest position and true 360 degree turning radius so you can keep a low center of gravity for moving heavy patients safely.

Additional accessories facilitate smooth transfer and transport of the patient, including the **Manta Transfer Sheet, the Bariatric Tow Package, and Push-Pull Handles** (see accessory page)

COMPLIANCE

SAE J3027 WHEN COMBINED WITH CORRESPONDING COT FASTENING SYSTEM


FERNO reserves the right to change specifications without notice.

© 2025 FERNO | WWW.FERNO.COM | U.S. 1-877-733-0911 | INTL. 1-937-382-1451 | NOV-1425-V5.1

BARIATRIC ACCESSORIES

LBS JR. BARIATRIC BOARD KIT

NO. 000LBSJRINX (INX)
NO. 000LBSJRPX1 (POWER X1)
NO. 000LBSJRPX2 (POWER X2)
NO. 0822184 (35X)



SX SURFACE XTENDER™

NO. FWESX (INX)
NO. 1590159 (POWER X1)
NO. 0822547 (POWER X2)



LBS JR. MATTRESS

NO. 0314050



LBS JR. CARRYING CASE

NO. 0822179 (MATT ONLY)
NO. 0822180 (MATT & LBS JR)



PILLOW WEDGE

NO. 0832159




MANTA TRANSFER SHEET

NO. 0382969



BARIATRIC TOW RING

NO. 0822511 (POWER X1)



BARIATRIC TOW RING

NO. 0822121 (35X)



IV POLE (LBS JR.)

NO. 0087173



STRAIGHT PUSH/PULL HANDLE

NO. 0822148 9HEAD END PX1 & 35X)
NO 0822112 (FOOT END PX1 & 35X)



ANGLED PUSH/PULL HANDLE

NO. 0822147 HEAD END PX1 & 35X)
NO 0822112 (FOOT END PX1 & 35X)



© 2025 FERNO | WWW.FERNO.COM | U.S. 1-877-733-0911 | INTL. 1-937-382-1451 | NOV-0425-V5