PEDIATRIC
TRANSPORT SOLUTIONS



Designed for the smallest patients.

Patients come in all shapes and sizes, and our smallest patients create unique challenges for patient transport. The PediMate®+, NeoMate™, and KangooFix™ provide efficient ways to safely secure pediatric patients from 3.5 to 100 lbs. to almost any cot without having to purchase additional equipment or seats.

PEDIMATE +
10 lb - 100 lb

NEOMATE
5 lb - 14 lb

KANGOOFIX
3.5 lb - 11.1 lb

FERNO

PEDIATRIC TRANSPORT SOLUTIONS



PEDIMATE +

- Fully-adjustable, five-point harness system securely holds patients, providing safe restraint for transport
- Rolls compactly for convenient storage and easy deployment
- Vinyl construction is non-toxic and easy-to-clean
- Three restraint straps easily attach to any cot

The PediMate+ provides a quick and efficient way to safely secure young patients to almost any cot without having to purchase additional equipment or specialized seats. Fully adjustable with a five-point harness system that securely holds patients, providing comfortable restraint for transport.

SPECIFICATIONS

ITEM NUMBER*
0314115

WEIGHT
3.5 lb / 1.59 kg

WIDTH
19" / 48.26 cm

LENGTH
35.5" / 90.17 cm

PATIENT WEIGHT
10 lb - 100 lb (4.5 kg - 45.4 kg)





NEOMATE

- Halo Pad™ keeps infant's head stable
- Fully-adjustable, five-point harness system fits children firmly and provides safe restraint during transport
- Three restraint straps easily attach to any cot
- Rolls up compactly for efficient storage and deployment

The NeoMate adapts an ambulance cot to safely transport pediatric patients and is fully adjustable with a five-point harness system that securely holds patients, providing comfortable restraint for transport.

SPECIFICATIONS

ITEM NUMBER*
0314116

WEIGHT
2.6 lb / 1.18 kg

WIDTH
10.5" / 26.67 cm

LENGTH
21" / 53.34 cm

PATIENT WEIGHT
5 lb - 14 lb (2.3 kg - 6.4 kg)





KANGOOFIX

- EN 1789 compliant for safety
- Skull cap included to help maintain baby's body heat
- Quick release system allows fast access for EMS professionals
- "Face-up" mode enables EMS professionals to quickly monitor newborn

The KangooFix Neonatal Restraint System was developed to safely and effectively cradle and secure a newborn during non-critical ambulance transportation, allowing the baby to travel with their parent.

SPECIFICATIONS

ITEM NUMBER*
0374906

WEIGHT
2.2 lb / 1 kg

WIDTH (MAX)
36.6" / 93 cm

LENGTH (MAX)
43.3" / 110 cm

PATIENT WEIGHT
3.5 lb - 11.1 lb (1.6 kg - 5 kg)



NASEMSO pediatric transport recommendations call for all devices, or combination of devices, to cover a weight range of between 5 to 99 pounds (2.3 - 45 kg). Combining the KangooFix, NeoMate, and PediMate + creates a covered weight range of 3.5-100 lb, supporting the safest transport possible for a wide range of ages and sizes.